

FlexiFuge - a new spin on things!



This is a versatile, personal micro-centrifuge for quick, low speed centrifugation.

Whisper quiet motor offers two speeds (2000g or 5000g / 6,000 or 9,200 RPM) depending on application.

FlexiFuge can be operated in continuous mode, pulsed for a short burst of rapid centrifugation, or timed to spin from 1-15 minutes. Regardless of setting, rotors stop once the lid release button is pushed.



LED lights as optical control illuminates when the unit is active and glow Red, Green or Blue. The color can be changed or turned off by the switch on bottom of unit.

FlexiFuge is supplied with 3 rotor's:
One 8-place fixed angle rotor for 1.5ml tubes and 8 adapters for 0.5ml tubes, one rotor for (4) 8 x0.2ml PCR tube strips and one 4-place rotor to spin 5ml micro-centrifuge tubes.



The features you've come to expect:



- Whisper quiet motor
- Two speeds - 6,000 or 9,200 RPM's
- Operates in continuous mode, pulse or timer (1-15 min.)
- NEW: rotor for 5ml !



FlexiFuge	Stock Item / Swiss warehouse	Unit Qty
AG.C1000E	FlexiFuge complete with 3 rotors; 1 rotor for 1.5 tubes + 8 adapters for 0.5ml tubes, 2 rotors 0.2ml PCR-strips & 5ml tubes	1 pk