

RotoFlex & RotoFlex Plus

New Item!

The **RotoFlex** will accommodate a variety of tubes from 1.5 – 50ml by simply changing the rotisserie. No tools are required! The mixing angle can be adjusted from 0-90° and each of the 4 tube panel can be aligned independently. The 360° rotation ensures thorough mixing of all samples at the same time while rotating at 22 rpm.

To change the rotation direction – simply touch the rotisserie paddle while in motion and the RotoFlex automatically begins to rotate in the opposite direction.

RotoFlex Plus has the additional feature of digital display, variable speed from 10-40rpm and oscillating mode at 40° rotation, controlled by an easy keypad. Optional there is a magnetic platform to optimize Western Blot hybridization. Simply attach your plastic bag, containing membrane & antibody solution with the metal strips provided.



RotoFlex and RotoFlex Plus can be used in incubators and cold room settings. The unit is supplied complete with base unit and rotisseries (1 pack of each size).

Ideal for:

- Hybridization
- DNA Purification
- Sample mixing / re-suspension
- Protein binding
- Immunoprecipitation



Digital face panel, with variable speed keypad (RotoFlex Plus model)



Removable Inserts - easily snap in and out of place



Change inserts to accommodate varying tube sizes



Simply touch to engage auto reverse feature



RotoFlex	Stock Item / Swiss warehouse	Unit Qty
AG.R2001	RotoFlex 220V, complete system	1 ea
AG.R2201	RotoFlex Plus 220V	1 ea
AG.R2012	Magnetic platform for RotoFlex Plus, Metal strips incl. (4)	1 ea
AG.R2013	Metal strips for magnetic plate for RotoFlex Plus, 4/pk	4 pk
AG.R2002	Rotisserie for 0.5 - 0.8ml tubes (21 tubes)	2 pk
AG.R2003	Rotisserie for 1.5 - 2ml tubes (15 tubes)	2 pk
AG.R2004	Rotisserie for 5 - 7ml tubes (7 tubes)	2 pk
AG.R2005	Rotisserie for 10 - 15ml tubes (7 tubes)	2 pk
AG.R2006	Rotisserie for 50ml tubes (6 tubes)	2 pk