

Omega PLUS™ Pipettor

Lightweight, ergonomic, UV resistant serological pipettor for accurate reproducible filling and dispensing of glass and plastic pipettes from 1 -100mL.

The **Omega PLUS™** Pipettor is a lightweight (190g), ergonomically designed pipette controller that combines single-handed pipetting with fast, precise liquid dispensing. The aspirating and dispensing speeds can be continuously adjusted with one hand using the thumb dial. Touch sensitive contoured buttons enable precise filling and dispensing with one hand. Its compact design makes it ideal for pipetting liquids in confined spaces.



Rechargeable 6.5V lithium ion battery will provide 10 hours of continuous use. **NEW: can be opened to replace old battery** by the user. Unique five bar LED display indicates battery strength and when to recharge. Overfilling and contamination are prevented by a 0.45µm hydrophobic filter located in the autoclavable nosecone. The Omega will accept glass or plastic pipettes from 1 – 100mL and is available in six colours. Unit is supplied with 6.5V lithium ion battery, 220V, 50Hz adapter, spare 0.45µm filter, wall mount, utility knife and spare I.D. labels.

Item No #	Omega Pipettor Stock Item / Swiss warehouse	Unit Qty
AG.P5015	Omega PLUS Green	1/EA
AG.P5016	Omega PLUS Red	1/EA
AG.P5017	Omega PLUS Blue	1/EA
AG.P5018	Omega PLUS Clear	1/EA
AG.P5019	Omega PLUS Purple	1/EA
Accessories		Unit Qty
AG.P1020	Silicone pipettor mount for pipettes 0.1 ml - 50 ml	1/EA
AG.153015	Filter, sterile, pore size 0.45 µm	10/PK
AG.153016	Filter, non-sterile, pore size 0.45 µm	10/PK
AG.P3225	Rubber filter interface	1/EA
AG.P5025	Nosepiece (Red)	1/EA
AG.P5070	Nosepiece (Grey)	1/EA
AG.P3010	Battery Charger	1/EA
AG.P1005	Pipettor Wall Mount	1/EA
AG.P5022	Omega I.D. Labels	5/PK