

e-Loop Plant

Annealing without flame

Ideal for instruments:
like Forceps and scissors



No.	Text E
LF5.101.000	e-Loop Large, 230V, NEW Diameter 21mm, length 161mm, with 7 angular adjustments, Residual heat display, Tray for 3 instruments, Screwdriver, Instruction manual and 2-Year Warranty
LF5.102.000	e-Loop Large, 120V, NEW Diameter 21mm, length 161mm, with 7 angular adjustments, Residual heat display, Tray for 3 instruments, Screwdriver, Instruction manual and 2-Year Warranty
LF5.103.000	e-Loop Large, 100V, NEW Diameter 21mm, length 161mm, with 7 angular adjustments, Residual heat display, Tray for 3 instruments, Screwdriver, Instruction manual and 2-Year Warranty
LF5.103.102	e-Loop stainless steel ring-holder for scalpel / forceps/scissors
LF5.103.200	e-Loop Timer + Holder, Stainless Steel

